15 mm 138 mm 3 mm

10 L
Purisoletm SM diluted

UK - Purisoletm SM diluted, Solution for Irrigation
1 L = 27,0 g Sorbitol, 5,4 g Mannitol
Air-vent needle behind booklet.

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Seite 12

UK

Purisoletm SM diluted should be used with caution in patients with diabetes mellitus (perioperative blood glucose monitoring recommended) or metabolic disorders (measuring of lactate levels recommended in the perioperative period). Adequate precaution measures should be taken to avoid massive fluid absorption, e.g. balancing of the fluid in- and outflow (stopping the procedure if the difference exceeds 1 - 2 liters), limiting the operation time to 1 hour, and not exceeding the recommended maximum pressure for the intervention.

Side-effects
Fluid absorption may occur in operative transurethral, hysteroscopic, and arthroscopic procedures. In case of intravascular absorption of large amounts of electrolyte-free irrigation solutions, hyperhydration and changes in the electrolyte balance (hyponatremia) may occur and consequently lead to acid base imbalance (metabolic acidosis). This condition is known as "TUR syndrome". The associated symptoms include nausea, vomiting, increase or decrease in blood pressure, bradycardia, headache, restlessness, confusion, lethargy, chest pain, muscular twitches, dyspnea, and poor urine output. In severe cases seizures, shock, pulmonary edema, cerebral edema, cardiovascular collapse, and coma may occur. Fatal cases have been reported. If the development of TUR syndrome is suspected, the procedure should be terminated, and adequate therapeutic measures initiated according to the individual condition of the patient. Potential side effects associated with the metabolism of sorbitol are hyperglycemia in patients with diabetes mellitus and hyperlactatemia and/or lactic acidosis in patients with metabolic disorders or in case of massive absorption.

Interaction with other medicines
None known.

Dosage
Dosage and quantity to be applied will result from the type of intervention.

Seite 3

UK INSTRUCTION FOR USE - READ CAREFULLY!
Purisoletm SM diluted - Solution for Irrigation Fresenius Kabi AG
61346 Bad Homburg
Germany

Composition
1 litre contains:
Sorbitol 27,0 g
Mannitol 5,4 g
Water for injections in bulk, Ph. Eur.
pH value 4,5 - 7,0
Titration acidity < 1 mmol NaOH/l
Theoretical osmolality 178 mosm/l

Contraindications
- Hypersensitivity to sorbitol and/or mannitol
- Hereditary fructose intolerance: Purisoletm SM diluted contains sorbitol. For this reason, the presence of hereditary fructose intolerance or a lack of fructose 1-6 diphosphatase must be ruled out prior to irrigation with Purisoletm SM diluted. Should anamnesistic data be unavailable or unclear, a fructose intolerance test must be performed before irrigation. Non-observation of this contraindication may lead to potentially fatal complications.
- Hyperhydration
- Methanol intoxication
- Newborns and infants

Precautions
Caution is required in patients with pre-existing hyponatremia. Serum sodium levels must be corrected before use in this case. Purisoletm SM diluted should be used with caution in patients with cardiac or renal insufficiency.

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Seite 5

UK

Descriptions and connectivity
The Kabi Container is a ten-liter canister. Cap and stopper can either be completely removed for pouring the solution or an irrigation transfer set can be spiked into the stopper of the Kabi Container. Based on tests conducted the following Fresenius Kabi irrigation transfer sets are compatible with the Kabi Container: Bladder Syringes Filling Set (Ref 2841001), Transfer Set 1D (Ref 2883101), Transfer Set 2 D/X (2877201), Twin Transfer Set 7 D (2890501), AT1 Auto-transfusion Set (9005101), Cystoscopy Transfer Set (Ref 2883131).

Use only if the solution is clear and free from particles and if the package is undamaged. For single use only. Any unused portion of the solution must be discarded. Re-use of medical single use equipment involves a risk of contamination of the equipment. This may cause infection in the patient or the user, which may lead to illness or death. Keep out of the sight and reach of children! Do not use Purisoletm SM diluted after the expiry date.

Date of revision
May 2021

Presentation
Polypropylene container (Kabi Container)
1 x 10,000 ml

Additional information
Not for intravenous use or injection!
The sterility of Purisoletm SM diluted is guaranteed only until the container is opened. The solution must be used immediately after opening.

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